

Compact Scales

HIL

Resolution **0.01g**

Low profile & small footprint

High Resolution

(1/10,000:HL-100, 1/4,000:HL-400/4000, 1/2,000:HL-200/2000)

Up to 1/10,000 resolution

HL-100

100g x 0.01g
(3.5oz x 0.0005oz / 500ct x 0.05ct)

HL-200

200g x 0.1g

HL-400

400g x 0.1g

HL-2000


2000g x 1g

HL-4000

4000g x 1g



HL-100



Standard carrying case

(HL-100)



HL-200



HL-2000

- For wide range of applications including schools, jewelry shops, offices and laboratories
- Low profile & small footprint - ideal for use in areas with space problems
- Portable operation using 6 AA batteries (not included) and AC adaptor (optional)
- Auto Power Off to preserve battery life
- High accurate weighing with adjustable feet and spirit levels
- HL-100 weighs in units of g. ct. and oz. (oz. only available in non-metric countries)
- Separate Power and Re-zero keys to eliminate operating errors

Stainless steel pan for sanitary weighing

(HL-100: Standard, HL-400/2000/4000 Optional)

| Model | HL-100 | HL-200 | HL-400 | HL-2000 | HL-4000 |
|------------|-----------------------|-------------------------------------|--------|---------|---------|
| Capacity | 100g/3.5oz/500ct | 200g | 400g | 2000g | 4000g |
| Resolution | 0.01g/0.0005oz/0.05ct | 0.1g | 0.1g | 1g | 1g |
| Pan Size | Ø70mm / Ø2.8inches | 136(W) x 136(D)mm / 5.4 x 5.4inches | | | |

Option

OP-10-HL Stainless Steel Weighing Pan
(For HL-400/2000/4000)

Accessories

AX-TB202 AC Adapter AC120V 'A' plug
AX-TB125 AC Adapter AC200-240V 'C' plug
AX-TB126 AC Adapter AC200-240V 'BF/S/NO' plug

ALYS Labware, Lausanne

Tel : 021.312.42.60

Fax : 021 312 42 61

labware@alys-technologies.com