

CMB Baby Scales



“The Affordable Alternative”

Adam Equipment Company Ltd
Bond Avenue
Milton Keynes, MK1 1SW
United Kingdom
Tel: +44 (0)1908 274545
Fax: +44 (0)1908 641339
e-mail: sales@adamequipment.co.uk

Adam Equipment Inc.
26 Commerce Drive
Danbury, CT 06810
United States of America
Tel: +1 203 790 4774
Fax: +1 203 792 3406
e-mail: sales@adamequipment.com


Adam Equipment S.A. (Pty) Ltd
P. O. Box 1422, Kempton Park 1620
Johannesburg
Republic of South Africa
Tel: +27 (0)11 974 9745
Fax: +27 (0)11 392 2587
e-mail: sales@adamequipment.co.za

CMB Baby Scales

“The Affordable Alternative”

The CMB 20 baby scale is a lightweight, portable scale for weighing both babies and toddlers.

- Easy to read LCD display in kg and lbs
- Water resistant keyboard and construction
- 9V battery or AC adapter powered (adapter supplied)
- Auto power off
- Hold function
- Hygienic and easy to clean
- Indicator can be wall mounted
- Tare function



SPECIFICATIONS

Model	CMB 20
Maximum capacity	20kg
Minimum capacity	200g
Readability	10g
Tare range	Full
Repeatability	+/- 1d
Linearity (±)	+/- 2d
Units of measure	kg, lb
Operating temperature	5 - 35°C
Power supply	9V dry cell battery or 9VDC adapter
Calibration	Automatic external
Display	1.0 LCD with 4 active digits
Pan size	520x265mm / 20.5x10.4"
Gross weight	2.5kg / 5.5lb
Applications	Baby / toddler weighing
Functions	Weighing
Other features	Removable cradle so that the unit can weigh small children, auto power off, display hold function, flexible cable to indicator with wall mount feature

Distributed by: [ALYS Labware, Lausanne](#)

Tel : 021 312 42 60

Fax : 021 312 42 61

labware@alys-technologies.com



Offices in UK, USA and South Africa

www.adamequipment.com