

MRP Personal Scale



“The Affordable Alternative”

Adam Equipment Company Ltd
Bond Avenue
Milton Keynes, MK1 1SW
United Kingdom
Tel: +44 (0)1908 274545
Fax: +44 (0)1908 641339
e-mail: sales@adamequipment.co.uk

Adam Equipment Inc.
26 Commerce Drive
Danbury, CT 06810
United States of America
Tel: +1 203 790 4774
Fax: +1 203 792 3406
e-mail: sales@adamequipment.com

Adam Equipment S.A. (Pty) Ltd
P. O. Box 1422, Kempton Park 1620
Johannesburg
Republic of South Africa
Tel: +27 (0)11 974 9745
Fax: +27 (0)11 392 2587
e-mail: sales@adamequipment.co.za

MRP Personal Scales


“The Affordable Alternative”

The MRP 200P personal scale is a fast, accurate scale for person weighing. The MRP 200H has the same features as the MRP 200P, but has a two metre height measuring pole in addition.

- Simple to use
- Platform has non-slip rubber surface
- Powered by internal batteries or external adapter (adapter supplied)
- Hold function to lock the display
- Auto off to maximise battery life
- Stability filters for person weighing
- Tare function

SPECIFICATIONS

Model	MRP 200P	MRP 200H
Maximum capacity	200kg	200kg
Minimum capacity	2kg	2kg
Readability	100g	100g
Resolution	100g	100g
Tare range	Full	Full
Repeatability	100g	100g
Linearity (±)	200g	200g
Units of measure	kg, lb	
Operating temperature	0°C - 40°C	
Power supply	6 alkaline batteries, AC adapter	
Calibration	Automatic External	
Display	LCD	
Pan size	250x300mm / 9.8x11.8"	
Overall dimensions (wxdxh)	320x400x1100mm / 12.6 x15.7x43.3"	
Note:	Height of MRP 200H is 2120mm / 83.5" due to height measuring pole	
Applications	Person Weighing	
Functions	Weighing	
Other features	Hold stability function, 2m height pole on MRP 200H	



Distributed by: [ALYS Labware, Lausanne](#)

Tel : 021 312 42 60
 Fax : 021 312 42 61
labware@alys-technologies.com



Offices in UK, USA and South Africa

www.adamequipment.com