



Sky Line

MAGNETIC STIRRER MS-01 AND VORTEX MIXER V-3

MAGNETIC STIRRER MS-01

MAGNETIC STIRRER MS01 is a powerful and simple to use four position stainless still magnetic stirrer.

PRODUCT FEATURES

- Four positions for effective mixing in different size flat bottom laboratory flasks.
- Compatibility with ELMI thermostat TW-2.03
- Stainless still highly resistible body.
- Can be easily integrated in different lab processes.

PRODUCT BENEFITS

It is specially designed to operate with ELMI water thermostat TW-2 series. Together with TW-2.03 heating functions it becomes possible to maintain highly accurate temperature maintenance as well as mixing efficiency control. Reliability, simple operation and noiseless operation make MS01 perfect assistant in any laboratory.



MAGNETIC STIRRER MS01 SPECIFICATIONS

Speed range	RPM	80 - 1 500
Liquid volume range for one mixing place	L	from 0.2 up to 2L
Power adapter output	⊖ ⊕	DC in 12V, 1 000 mA
Dimensions (length x width x height)	mm	210 x 270 x 35
Weight	kg	1.8

VORTEX MIXER V-3

VORTEX V-3 is a pretty looking and very easy handled vortex mixer for any type of test tubes of up to 50 ml volume.

PRODUCT FEATURES


- Touch or hands free mode.
- Easy vortex intensity regulation.
- 2 exchangeable vortexing adapters.
- Quiet and maintenance free mechanism.
- Half-sphere shape for maximal stability on the desktop.

PRODUCT BENEFITS



2 TYPES OF EXCHANGEABLE VORTEX ADAPTERS

CONICAL-TYPE ADAPTER GRID-TYPE ADAPTER



ALYS Labware, Lausanne

Tel : 021 312 42 60

Fax : 021 312 42 61

labware@alys-technologies.com

VORTEX MIXER V-3 SPECIFICATIONS

Socket oscillation speed	RPM	100 to 2500
Socket oscillation amplitude	mm	3
Power adapter output	⊖ ⊕	DC in 12V, 500 mA
Power consumption	W	5
Dimensions (diameter x height)	mm	130 x 60
Weight	kg	0.6